

# Easyboot Epic & Grip Application



1.

**1. PUTTING ON THE EASYBOOT EPIC**  
Untie the 3/4" nylon webbing pull-strap\*\* that is tied around the heelstrap. Loop it around the heelstrap so that it is easily removed once the boot is applied. Fold the gaiter over the back of the boot before fitting the boot over the toe. Loosen the cable and make sure the buckle is open. Pick up the horse's leg and support it between your knees so that you can use two hands when applying the boot. Place the boot over the horse's toe and, using a slight back and forth twisting motion, slide the boot on the rest of the way keeping the gaiter from bunching into the boot.



1b.



2.

**2. USING THE PULL-STRAP\*\***  
After the boot is on as far as it will go, use the nylon pull-strap (one included with each boot) for leverage to pull the boot the rest of the way on, as well as to pull the inside heelstrap out from under the sole of the hoof. Make sure the heelstrap is not pulled into the soft area of the heel bulbs. If you are applying the boot over steel shoes, make sure that the heelstrap is above the back of the shoe.



4.

\*\*In many cases you will not need to use the nylon pull strap. If you do not use it, be careful not to pull on the gaiters. Make sure you remove the nylon pull strap completely from the boot prior to riding.




4.

**3. SETTING THE HEEL**  
Put the foot down on the ground and allow the horse to put weight into it. From this position, you may have better leverage with the pull-strap. Remove the pull-strap after the boot has been successfully applied.

**IMPORTANT!**  
Please read this section carefully before fastening your gaiter.


**4. FASTENING THE EASYBOOT GAITER**  
Pull the top of the gaiter tight around the pastern. Fasten the hook and loop straps. The tighter the fit, the better the boot will perform.

Due to their great deal of elasticity, gaiter straps cannot be over tightened. Overlap the Velcro as much as needed to ensure a tight fit and avoid rubbing issues. This overlapping does not compromise the integrity of the gaiter.



When properly tightened, it should be difficult to fit a finger between the pastern and the top of the gaiter.

**5. LOCKING THE CABLE AND BUCKLE**  
After the boot is on, the inside heelstrap is in the correct position, and the gaiter is fastened, put the horse's foot down and set the cable in the tightest position possible. Clamp the buckle down either (a) using the palm of your hand or (b) using your foot. When tight and closed correctly, the buckle should be difficult to open without the aid of a screwdriver or hoof pick. Although the buckle has three adjustments, additional cable tightness can be achieved with different cable wraps (see diagrams). Make sure to unclamp the buckle before adjusting cable.



- a. **Standard cable:** This is the normal position of the cable at the top of the buckle.
- b. **Loose cable wrap:** If more slack is needed, adjust the cable into this position.
- c. **Tight cable wrap:** If the cable needs to be tightened, position the cable in this loop formation, around the top of the buckle.
- d. **Incorrect cable wrap:** Wrapping the cable in this way will result in frayed or broken wires. The buckle will also be more likely to pop up during use. Do NOT wrap the cable around the buckle clamp.

**6. INSERT THE OPTIONAL EASYBOOT BUCKLE PIN**  
The Easyboot Buckle Pin can now be added to lock the buckle into place. The Easyboot Buckle Pin is recommended for rocky trails, endurance events or on horses that have a tendency to open the front Easyboot buckle.

**7. A PROPERLY APPLIED EASYBOOT EPIC**  
The heelstrap is designed to fit snugly, low around the quarters, below the hairline, and below the heel bulbs. If the horse has a very low heel, make sure the strap is not pulled out of its position, or in any way allowed to abrade a low heel. Horses with low heels may require the heelstrap to be cut down a notch or two to avoid chafing near the heel bulbs. Do not use the boot if any portion of the heelstrap is above the hairline or touching sensitive skin.

**IMPORTANT:** The Easyboot Epic should fit tight around the hoof and the Gaiter snugly around the pastern. The Easyboot Epics should require effort to put on. After you fit the Easyboot Epic and the horse has had a chance to move around, make sure the Easyboot Epics are still fitted squarely and securely on the hoof. Make sure the Easyboot Epic has not twisted to the side. If the Easyboot Epic has twisted to the side it is either too large or needs to be tightened. Discontinue use if the boot continues to twist. Boots that twist can cause damage to the hoof wall.

For flexibility in cold weather, keep the Easyboot Epic in a heated room or soak in hot water before fitting on hoof.

Normal use of Easyboot Epics implies removal when not required. It is NOT recommended that the Epic be left on the hoof for more than 24 hours.

## 8. REMOVING THE EASYBOOT EPIC OR GRIP

The easiest way to remove a properly fitted Easyboot is to use a flat-head screwdriver. (a) Open the buckle with the screwdriver or hoofpick and unhook the Gaiter straps. (b) Next, place the screwdriver inside the boot, between the boot and the hoof at the quarter. Wedge the screwdriver under the hoof to pry the boot off. You may have to work both sides of the boot for easier removal. If the boot is fitted over steel shoes, you will have to work the heelstrap over the back of the shoe for removal.



6.



7.



8a.



8a.



8b.



# EasyCare Inc. Warranty

EasyCare offers a 90 day repair or replacement warranty for all hoof boots (purchased as new) manufactured under the EasyCare label. Please register your product at [http://www.easycareinc.com/warranty\\_reg.aspx](http://www.easycareinc.com/warranty_reg.aspx) and keep your original dated receipt(s) as proof of purchase during the warranty period.

If a boot fails during warranty, it is generally a minor hardware issue. Please bear in mind that most components of EasyCare hoof boots are either replaceable or repairable. This includes wire laces, cables, buckles, gears, rollers, screws and straps. Most repairs can be easily accomplished with simple household tools. Depending on the circumstances, EasyCare will send you repair parts, or ask that you return the boots to us for evaluation. Videos of many repairs are available on our website at [www.easycareinc.com](http://www.easycareinc.com).

While this warranty is very broad in nature, it does not cover wear (including sole breakthrough when worn over horseshoes), misuse, abuse, factory seconds, bargain bin items, or boots purchased as used. Products returned for evaluation must be clean, dry and free of debris, dirt, sand and manure. *Failure to return clean products will result in product being shipped back at sender's expense or a cleaning fee may be assessed.* Damage inflicted by omnivorous horses, rodents or other unusual abuse is not covered.

Returns must be authorized before sending and must be returned with an RMA number. Contact EasyCare directly for all warranty issues at 1-800-447-8836 to receive an RMA number.

Specifically:

1. Returns for credit only (refund of purchase price only): Products must be in new, unused condition, free of dirt and wear marks. It is your responsibility to evaluate the fit of your new product in a manner where you may return it in new condition if you have doubts about the fit. Be certain that you know the return policy of product purchased from dealers other than direct from EasyCare (the manufacturer). EasyCare will not give credit for product purchased elsewhere. Return to point of purchase in that circumstance. Also, all brochures, inserts and/or free gifts must accompany returned items.

Online Warranty Registration:  
[www.easycareinc.com/warranty\\_reg.aspx](http://www.easycareinc.com/warranty_reg.aspx)

2. Exchanges of new product (size or color issue): EasyCare will only exchange NEW items with same product style and model for same product when requested with NEW items only. You may exchange for a different size or color exchange only, no exceptions. If you wish to have a different style of hoof boot, etc., you must return new item for refund (original point of purchase) and then purchase new choice.

3. Warranty failure (used product that is within 90 days of purchase): Contact us at 800-447-8836. EasyCare will either provide repair parts directly to you, or repair/replace after our evaluation. Items must be returned to EasyCare at the consumers cost with original receipt. EasyCare will return the items via ground service (UPS or FEDX) free of charge. Requests for faster shipping service will be at the expense of the consumer.

4. Not under warranty (over 90 days or no proof of purchase): Repair parts may be purchased from your local dealer or EasyCare, or product may be returned for evaluation and possible repair (parts, labor and freight charges will apply).

5. Call for an RMA number (under warranty or not): Returns must be authorized before sending and must be returned with an RMA number clearly written on outside of packaging. Contact us at 800-447-8836 to receive an RMA number. Returns will be refused without an RMA number on outside of package.

6. Freight charges (outside warranty period): Freight charges to and from EasyCare for repairs, exchanges and items outside the warranty period will be the purchaser's responsibility.

7. EasyCare staff is always willing to discuss any specific issues you may have regarding the performance of your EasyCare or Old Mac's products. Please contact us at 800-447-8836 or email to [admin@easycareinc.com](mailto:admin@easycareinc.com). Visit our website, [www.easycareinc.com](http://www.easycareinc.com), for tips and product information that will help with the use of any of our product line.



**Easycareinc.com 800.447.8836**

Additional help is only a click away! Visit these links to learn more.

- 1. Find a dealer in your area**  
[www.easycareinc.com/Search/Dealer.aspx](http://www.easycareinc.com/Search/Dealer.aspx)
- 2. Find a Hoof Care Professional**  
[www.easycareinc.com/Search/Practitioner.aspx](http://www.easycareinc.com/Search/Practitioner.aspx)
- 3. Get expert hoof boot fitting advice and the inside scoop from our blogs**  
[www.easycareblogs.com](http://www.easycareblogs.com)
- 4. Read the latest in hoof boot and hoof care articles**  
[www.easycareinc.com/education/articles.aspx](http://www.easycareinc.com/education/articles.aspx)
- 5. Check out our instructional videos**  
[www.easycareinc.com/education/videos.aspx](http://www.easycareinc.com/education/videos.aspx)
- 6. Subscribe to our free monthly e-newsletter**  
[www.easycareinc.com/NewsRoom/Archives.aspx](http://www.easycareinc.com/NewsRoom/Archives.aspx)
- 7. Download product brochures, manuals and instructions**  
[www.easycareinc.com/education/downloads.aspx](http://www.easycareinc.com/education/downloads.aspx)
- 8. Check out our contests and promotions**  
[www.easycareinc.com/Cool\\_Staff/promotions.aspx](http://www.easycareinc.com/Cool_Staff/promotions.aspx)
- 9. Read our press releases for breaking news**  
[www.easycareinc.com/NewsRoom/Media.aspx](http://www.easycareinc.com/NewsRoom/Media.aspx)
- 10. Get discount prices at our Bargain Bin**  
[www.easycareinc.com/externallink.asp?C=XBAR](http://www.easycareinc.com/externallink.asp?C=XBAR)
- 11. Trade in boots at the Hoof Boot Upgrade Program**  
[www.hoofbootupgrades.com](http://www.hoofbootupgrades.com)
- 12. Team Easyboot**  
[www.teameasyboot.com](http://www.teameasyboot.com)
- 13. Learn how to do your own hoof boot repairs**  
[www.easycarerrepairs.com](http://www.easycarerrepairs.com)

We continue to strive for the best products on the market with you and your horse in mind. We hope you enjoy our product selection and website improvements. Thank you for your continued support and patronage.

#### The Fine Print:

EasyCare Inc. hereby disclaims responsibility for the misuse of any products described in this brochure and web site and/or losses incurred through the use of products by any method other than that which is outlined in this catalog or web site or EasyCare's product packaging. All products should be used in accordance with the instructions printed on the product's label. EasyCare Inc products are not intended to treat, cure or diagnose any medical condition.



**Easyboot® Epic**  
and  
**Easyboot® Grip**  
Application Guide

## Comfort Pads

Give your horse more comfort in each step! Comfort Pads are available in two thicknesses to allow you to decide what level of extra cushioning you give your horse. The pads are available in the original 12mm thickness, as well as the newer, thinner 6mm thickness. This thinner pad option may be used by itself or stacked in any desired combination for your custom application and is great for modifying boot fit.



## Dome Comfort Pads

These pads have proven to be very popular with a variety of users including trail riders, endurance riders and hoof care professionals. We recommend that you work with your veterinarian or hoof care practitioner when treating any type of hoof ailment.



I use these pads in Bares, Epics and Easyboots every time I ride. They are perfect when transitioning a horse from shoes to barefoot. My horses love them!

—Diane Conklin

I absolutely love the comfort pads. They hold up very good and they offer my horse the frog stimulation she should have. They can also help you to get your boots to fit a little better when you just need a little extra help with the fit.

—Lisa Grigaitis

My horses love these pads. When I put on their boots with these pads in they really lick their lips when I set their foot down. After a long ride using a new set of pads, it is obvious looking at the pads how much support they have given to the horse's foot, as well as everything else up the leg. Good job EasyCare! —Kerry Grear

Have Questions? Need Help? Visit [www.easycareinc.com](http://www.easycareinc.com)



**EasyCare**  
[www.easycareinc.com](http://www.easycareinc.com)